



GeneMATRIX Bio-Trace DNA Purification Kit

Kit for isolation of DNA from biological traces

Cat. no. E3510

Version 3.1

January, 2008

For laboratory use only.
Not for drug, household or other uses.

Note 1: The Bio-Trace Kit is designed for isolation of DNA from the traces of biological samples, such as: fresh or frozen blood, blood, saliva, semen spots, hair roots and stems, cigarette filters, chewing gum, fragments of tissues, urine, among others.

Note 2: Once the kit is unpacked, store components at room temperature, with the exception of Sol BT buffer and Proteinase K. Sol BT buffer should be kept at $2\pm 8^{\circ}\text{C}$ and Proteinase K at -20°C .

Note 3: Tissue lysates are very sticky. This can lead to slow lysate filtration through the resin. Therefore it is advisable to check, if lysate and washes passed completely through the resin.

Note 4: All solutions should be kept tightly closed to avoid evaporation and resulting components concentration changes.

Note 5: The kit does not contain 1M DTT, 96 % ethanol, xylene and PBS. To prepare sterile PBS, dissolve 8 g NaCl, 0.2 g KCl, 1.44 g Na_2HPO_4 and 0.24 g KH_2PO_4 in 800 ml H_2O . Adjust pH to 7.4 with HCl. Add H_2O to 1 liter.

Protocol

I. DNA binding spin-columns activation

1. Apply 40 μl of activation **Buffer BT** onto the DNA binding spin-column (do not spin) and keep it at room temperature till transferring lysate to the spin-column.

Note 1: Addition of Buffer BT onto the center of the resin enables complete wetting of membranes and maximal binding of DNA.

Note 2: The membrane activation should be done before starting isolation procedure.

II. Sample preparation

A. 1-100 μl of blood

1. Add to the 2 ml Eppendorf tube:
1-100 μl of fresh or frozen blood (if the blood volume is less than 100 μl , fill up with **PBS** to 100 μl), 250 μl **Lyse BT** buffer and 10 μl **Proteinase K**.
2. Mix by inverting the tube several times or by vortexing.
3. Follow the point 1. Part III. of the DNA isolation protocol.

B. Spots of blood, saliva, semen and other biological liquids

- 1.a. Cut off the fragment of spotted material (do not exceed 1 cm^2). Cut the material into small pieces. Place the pieces in the 2 ml Eppendorf tube.
- 1.b. Spots from solid surfaces should be scratched and poured into the 2 ml Eppendorf tube.
2. Add 350 μl **Lyse BT** buffer and 10 μl **Proteinase K**. In the case of semen add additionally 20 μl **1M DTT**.
3. Mix thoroughly by inverting the tube and incubate for 60 min at 56°C . Mix by inverting every 15 min.
4. Follow the point 1. Part III. of the DNA isolation protocol.

C. Hair roots and stems, nails

- 1.a. Cut off 0.5-1 cm fragments of hair including the roots or cut hair stems into 0.5-1 cm long pieces (if the sample does not contain hair roots).
 - 1.b. Cut nails into small pieces (less than 2 mm²).
 2. Place cut hair or nails (maximum sample size: 20 mg) in 2 ml Eppendorf tube. Add 350 µl **Lyse BT** buffer, 20 µl **1M DTT** and 20 µl **Proteinase K**.
 3. Mix by inverting the tube several times or by vortexing, then centrifuge to collect hair or nails at the bottom of the tube. Incubate hair roots or hair stems for at least 60 min at 56°C (until completely digested). Mix by inverting the tube every 30 min. In the case of nails overnight incubation is recommended.
- Note 1:** While isolating DNA from hair sample containing hair roots, make sure the roots are completely submerged in the lysis buffer. If necessary, submerge the parts of hair containing roots in the lysis buffer using a sterile tips.
4. Follow the point 1. Part III. of the DNA isolation protocol.

D. Cigarette filters, chewing gum

- 1.a. Remove outer paper from the filter. Cut the paper into small pieces and place the sample in the 2 ml Eppendorf tube.
- 1.b. Chewing gum (maximum sample size: 50 mg) cut into small pieces and place in the 2 ml Eppendorf tube.
2. Add 350 µl **Lyse BT** buffer and 10 µl **Proteinase K**.
3. Mix by inverting the tube several times or by vortexing and incubate for 60 min at 56°C. Mix by inverting every 15 min.
4. Follow the point 1. Part III. of the DNA isolation protocol.

E. Fragments of tissues

1. Cut the tissue fragment (maximum sample size: 10 mg) into small pieces and place the sample in the 2 ml Eppendorf tube. Suspend the sample thoroughly in 350 µl **Lyse BT** buffer. Add 20 µl **Proteinase K**.
2. Mix by inverting the tube several times or by vortexing and incubate for 3-6 h at 56°C, mix by inverting every 15 min.

Note 1: After incubation the tissue suspension should be completely digested and appear as transparent fluid. If not, continue the incubation.

3. Follow the point 1. Part III. of the DNA isolation protocol.

F. Paraffin-embedded tissues

1. Prepare a small section (up to 25 mg) from block of embedded tissue. Place the sample in 2 ml Eppendorf tube.
2. Add 1 ml **xylene**. Vortex vigorously.
3. Incubate at room temperature for 15 min.
4. Centrifuge for 3 min at 12000 rpm.
5. Remove supernatant by pipetting.
6. Add 1 ml **xylene** to the pellet, mix by vortexing.
7. Centrifuge for 3 min at 12000 rpm.
8. Remove supernatant by pipetting.
9. Add 1 ml **96 % ethanol** to the pellet. Mix by vortexing or inverting the tube.
10. Centrifuge for 3 min at 12000 rpm.
11. Remove supernatant by pipetting.
12. Repeat steps 9-11 once.
13. Incubate the open tube at 37°C until the ethanol has evaporated (app. 15 min).
14. Resuspend the tissue pellet in 350 µl **Lyse BT** buffer. Add 20 µl **Proteinase K**.
15. Follow the point 2. of **E. Fragments of tissues** protocol.

G. Formalin-fixed tissues

1. Wash tissue sample twice with **PBS** to remove fixative. Discard **PBS**.
2. Cut tissue fragment (up to 25 mg) into small pieces. Place the sample in 2 ml Eppendorf tube and add 350 µl **Lyse BT** buffer and 20 µl **Proteinase K**.
3. Follow the point 2. of **E. Fragments of tissues** protocol.

H. Urine

1. Add 2 ml of urine to the 2 ml Eppendorf tube.
2. Centrifuge urine in microcentrifuge for 2 min at 8000 rpm.
3. Carefully discard the supernatant without disturbing the pellet. Add to the pellet 350 μ l **Lyse BT** buffer, 20 μ l **1M DTT** and 10 μ l **Proteinase K**.
4. Vortex for 15 sec.
5. Incubate for 60 min at 56°C, mix by inverting the tube every 15 min.
6. Follow the point 1. Part III. of the DNA isolation protocol.

III. DNA isolation

1. Add 350 μ l of **Sol BT** buffer and mix thoroughly by several times inverting the tube.
2. Incubate for 10 min at 70°C.
3. In the case of biological material samples mentioned in sections:
A do not add of **96-100 % ethanol**
B, D, H add 180 μ l of **96-100 % ethanol**
C add 230 μ l of **96-100 % ethanol**
E, F, G add 350 μ l of **96-100 % ethanol**
4. Mix thoroughly by several times inverting the tube.
5. Centrifuge for 2 min at 14000 rpm.
- 6.a. Transfer 600 μ l of supernatant to the **DNA binding spin-column**, placed in the collection tube.
- 6.b. In the case of isolation of DNA from dried spots (blood, saliva or semen) on fabric-type material, transfer the supernatant to the **DNA binding spin-column**, placed in the collection tube (do not exceed the total volume of 600 μ l of the supernatant). The soaked material with the remainings of the supernatant left in the tube after lysis should be transferred to the **filtration spin-column** to recover the entire lysate. Centrifuge **filtration spin-column** for 2 min at 12000 rpm. Save the flow-through for the following isolation procedure.
7. Centrifuge **DNA binding spin-column** for 1 min at 12000 rpm.
8. Remove the spin-column, discard flow-through and place back spin-column in the collection tube.

9. Transfer the remainings of supernatant (point 6.a.) or the flow-through from the **filtration spin-column** (point 6.b.) to the **DNA binding spin-column**, placed in the collection tube. Repeat centrifugation for 2 min at 12000 rpm to pass completely the lysate through the resin.

Note 1: Continue centrifugation, if not all of the lysate passed through the column.

10. Take out **DNA binding spin-column**, discard flow-through and place back spin-column in the collection tube.
11. Add 500 µl of **Wash BTX1** buffer to the **DNA binding spin-column** and centrifuge for 1 min at 12000 rpm.
12. Take out **DNA binding spin-column**, discard flow-through and place back spin-column in the collection tube.
13. Add 500 µl of **Wash BTX2** buffer to the **DNA binding spin-column** and centrifuge for 2 min at 12000 rpm.
14. Place **DNA binding spin-column** in a new collection tube (1.5-2 ml) and add 50 µl of **Elution** buffer (10 mM Tris-HCl, pH 8.5) heated to 70°C to elute bound DNA.

Note 1: Addition of eluting buffer directly onto the center of the resin improves DNA yield. To avoid transferring traces of DNA between the spin-columns do not touch the spin-column walls with the micropipette.

Note 2: The following eluting solutions can be used:

1. 5-10 mM Tris-HCl buffer, pH 8.0-9.0
2. 0.5-1 x TE buffer, pH 8.0-9.0 (not recommended for DNA sequencing).
3. Other special application buffers can be used, provided that their pH and salt concentration is similar to that of 5-10 mM Tris-HCl, pH 8.0-9.0.

15. Incubate **DNA binding spin-column/collection tube** assembly for 5 min at room temperature.
16. Centrifuge the **DNA binding spin-column** for 1 min at 12000 rpm.
17. Discard **DNA binding spin-column**, cap the collection tube. DNA is ready for analysis/manipulation. It can be stored either at 2÷8°C or at -20°C.

SELECTION OF THE KITS DEPENDING ON THE TYPE OF ISOLATED MATERIAL			AGAROSE - OUT	BACTERIAL & YEAST GENOMIC	BIO - TRACE	BASIC	BONE	CELL CULTURE	FOOD	HUMAN BLOOD RNA	PCR / DNA CLEAN-UP	PLANT & FUNGI	PLASMID MINIPREP	QUICK BLOOD	SHORT / DNA CLEAN-UP	SOIL	STOOL	SWAB EXTRACT	TISSUE	TISSUE & BACTERIAL	UNIVERSAL DNA/RNA /PROTEIN	UNIVERSAL RNA	UNIVERSAL RNA /miRNA	MICELLULA DNA
DNA	GENOMIC	BACTERIA		X																X				
		YEAST		X																				
		CELL CULTURE						X											X	X				
		PLANT AND FUNGI										X												
		BLOOD												X										
		SOIL														X								
		STOOL															X							
		SWAB																X						
		SOLID TISSUES																	X	X				
		LIQUID TISSUES																	X	X				
		RODENT TAILS																	X	X				
		HAIR																	X	X				
		INSECTS																	X	X				
		URINE																	X	X				
		BONE					X																	
		BIOLOGICAL TRACES			X																			
		FOOD							X															
	PLASMID	BACTERIA				X							X											
		YEAST		X																				
	ISOLATION FROM AGAROSE GELS		X			X																		
	PURIFICATION OF PCR PRODUCTS/DNA AFTER ENZYMATIC REACTIONS					X					X				X									X
DNA/RNA/PROTEIN FROM THE SAME BIOLOGICAL SAMPLE	ANIMAL TISSUE																				X			
	PLANT TISSUE																				X			
	BACTERIA																				X			
	YEAST																				X			
	CELL CULTURE																				X			
RNA	TOTAL RNA LONGER THAN 200 BASES	ANIMAL TISSUE																				X		
		PLANT TISSUE																				X		
		BACTERIA																				X		
		YEAST																				X		
		CELL CULTURE																				X		
		HUMAN BLOOD							X															
	miRNA AND TOTAL RNA	ANIMAL TISSUE																					X	
		PLANT TISSUE																					X	
		CELL CULTURE																					X	

GeneMATRIX is synthetic, new generation DNA- and RNA-binding membrane, selectively binding nucleic acids to composite silica structures. Novel binding and washing buffers are developed to take full advantage of **GeneMATRIX** capacity, yielding biologically active, high-quality nucleic acids. Matrix is conveniently pre-packed in ready-to-use spin-format. Unique chemical composition of the matrixes along with optimized construction of spin-columns improve the quality of final DNA or RNA preparation. To speed up and simplify isolation procedure, the key buffers are colour coded, which allows monitoring of complete solution mixing and makes purification procedure more reproducible.

As a result, we offer kits, containing matrixes and buffers that guarantee rapid, convenient, safe and efficient isolation of ultrapure nucleic acids. Such DNA or RNA can be directly used in subsequent molecular biology applications, such as: restriction digestion, dephosphorylation, kinasing, ligation, protein-DNA interaction studies, sequencing, blotting, *in vitro* translation, cDNA synthesis, hybridization among others. Additional advantage is reproducibility of matrix performance, as component preparation is carried at Eurx Ltd.

GeneMATRIX Bio-Trace DNA Purification Kit is designed for rapid purification of total DNA (genomic, mitochondrial) from the traces of biological samples (especially forensic case work samples), exemplified by: fresh or frozen blood, blood spot, saliva, semen spots, hair roots and stems, cigarette filters, chewing gum, fragments of tissues, urine. The kit is effective in DNA purification from both fresh samples and dried, many years old or preserved for example in alcohol or formaline. Purified DNA is free of contaminants, such as: proteins, lipids, dyes, detergents, organic inhibitors of enzymatic reactions, buffers, salts, divalent cations, among others. Biological trace sample is solubilized/lysed in the presence of special buffer, which desintegrates remaining tissue- and cellular structures, while preserving integrity and stimulating quantitative recovery of all traces of DNA spontaneously released during progressing decay of older samples. Further, Proteinase K digests contaminating proteins, including stripping-off DNA of all bound proteins, among them nucleases. Optimized buffer and ethanol are added to provide selective conditions for DNA binding during brief centrifugation, while contaminants pass through the **GeneMATRIX** resin in the spin-column. Traces of contaminants remaining on the resin are efficiently removed in two wash steps. High-quality cellular DNA is then eluted in low salt buffer, e.g.: Tris-HCl, TE or water. Isolated DNA is ready for downstream applications without the need for ethanol precipitation.