

Chill Rack qPCR (Soğutucu Rack)

- Özel akrilik, sızdırmaz materyalden üretilmiştir.
- 4 °C altında tahmini süre: 2,5 saat.
- Kullanımdan önce -20 °C' de 8 saat dondurulmalıdır.
- Otoklavlanamaz.

> NB-Chill01

- 64 x 0.1 ml Rotor tüp
- 6 x 0.5 ml tüp
- 6 x 1.5-2.0 ml tüp
- 96 x 0.2ml Plate
- Kapaklı kullanım




> NB-Chill02

- 4 x 1.5-2.0ml tüp
- 4 x 0.5 ml tüp
- 96 x 0.2ml Plate



> NB-Chill03

- 40 x 1.5-2.0ml tüp



+90 312 220 22 80

bm@bmlabosis.com

Kızılırmak Mah. 1441. Cad. Meva İş Merkezi No:41/8 A

www.bmlabosis.com

Chill Rack qPCR (Cooling Rack)

- Produced from special acrylic, leak-proof material.
- Estimated time under 4 °C: 2.5 hours.
- Must be frozen at -20 °C for 8 hours before use.
- Not autoclavable.

> NB-Chill01

- 64 x 0.1 ml rotor tube
- 6 x 0.5 ml tube
- 6 x 1.5-2.0 ml tube
- 96 x 0.2 ml plate
- Suitable for use with lid




> NB-Chill02

- 4 x 1.5-2.0 ml tube
- 4 x 0.5 ml tube
- 96 x 0.2 ml plate



> NB-Chill03

- 40 x 1.5-2.0ml tube



+90 312 220 22 80

bm@bmlabosis.com

Kızılırmak Mah. 1441. Cad. Meva İş Merkezi No:41/8 A

www.bmlabosis.com