

TipOne® Repeat Dispenser Tips

TipOne Repeat

Dispenser tips you'll want to use again!


TipOne® Repeat offers high quality, value replacement dispenser tips for many popular dispensing systems.

Features

- Available non-sterile bulk, or sterile individually wrapped
- High accuracy and precision
- ISO 8655 manufactured
- Polypropylene cylinder
- Polyethylene piston

TipOne Repeat Tips are compatible with:

- STARLAB RepeatOne® (see page 50)
- Eppendorf® Multipette® 4780
- Brand® Handystep® Manual
- Gilson® Distriman®
- Ritter® Ripette®
- Biohit® Electronic



Dispenser Tips		
Cat. No.	Description	Pack Size
Non-Sterile		
S4770-0050	0.5 ml TipOne® Repeat Dispenser Tips	100
S4772-0125	1.25 ml TipOne® Repeat Dispenser Tips	100
S4774-0250	2.5 ml TipOne® Repeat Dispenser Tips	100
S4776-0500	5.0 ml TipOne® Repeat Dispenser Tips	100
S4778-1250	12.5 ml TipOne® Repeat Dispenser Tips	100
S4780-2500	25 ml TipOne® Repeat Dispenser Tips	50
S4782-5000	50 ml TipOne® Repeat Dispenser Tips	25
S4026-1525	Adaptor for 25 ml and 50 ml Dispenser Tips	1
Sterile, Individually Wrapped		
S4784-0050	0.5 ml TipOne® Repeat Dispenser Tips (Sterile)	100
S4786-0125	1.25 ml TipOne® Repeat Dispenser Tips (Sterile)	100
S4788-0250	2.5 ml TipOne® Repeat Dispenser Tips (Sterile)	100
S4790-0500	5.0 ml TipOne® Repeat Dispenser Tips (Sterile)	100
S4792-1250	12.5 ml TipOne® Repeat Dispenser Tips (Sterile)	100
S4794-2500	25 ml TipOne® Repeat Dispenser Tips (Sterile)	50
S4796-5000	50 ml TipOne® Repeat Dispenser Tips (Sterile)	25

Sterile adaptors available on request.

Repeat Dispenser

See page 50
for ordering
details of
STARLAB's
RepeatOne®
Dispenser



Application Guide

Dispenser Tip Volume	Dispenser Adjustment				
	1	2	3	4	5
0.5 ml	10 µl	20 µl	30 µl	40 µl	50 µl
1.25 ml	25 µl	50 µl	75 µl	100 µl	125 µl
2.5 ml	50 µl	100 µl	150 µl	200 µl	250 µl
5.0 ml	100 µl	200 µl	300 µl	400 µl	500 µl
12.5 ml	250 µl	500 µl	750 µl	1000 µl	1250 µl
25 ml	500 µl	1000 µl	1500 µl	2000 µl	2500 µl
50 ml	1000 µl	2000 µl	3000 µl	4000 µl	5000 µl
Steps*	48	23	15	11	8

*Number of steps may vary according to dispenser manufacturer.