

PCR Sample Racks

Polypropylene provides high chemical resistance. Temperature use: -90 °C to 121 °C. Autoclavable.

PCR Workstation				
	Description	Colour	Cat. No	Pack Size
	PCR Workstation Polypropylene racks for the storage of PCR tubes, PCR strips and PCR plates, as well as 0.5 ml and 1.5/2.0 ml microcentrifuge tubes. Moulded, alpha-numeric matrix. The two different sized racks (Rack A and Rack B) can be combined to form the complete PCR Workstation. Racks can either be purchased separately or together as a complete workstation.			
	Rack A only (Satellite Rack) 96 wells to hold 0.2 ml PCR tubes, PCR strips or 1 x 96-well PCR plate Each rack supplied with transparent, lift-off lid Dimensions (WxDxH): 127 x 87 x 33 mm (with lid)	Neon Blue*	E2396-1641	1
		Neon Green*	E2396-1642	1
		Neon Orange*	E2396-1643	1
		Neon Pink*	E2396-1644	1
		Neon Yellow*	E2396-1646	1
		Mixed Neon**	E2396-9649	5
	Rack B only (Workstation Rack) Holds 32 x 0.2 ml PCR tubes, 16 x 0.5 ml tubes (PCR or MCT), and 24 x 1.5/2.0 ml microcentrifuge tubes Each rack supplied with transparent, lift-off lid Dimensions (WxDxH): 215 x 115 x 50 mm (with lid)	Neon Blue*	E2396-1651	1
		Neon Green*	E2396-1652	1
		Neon Orange*	E2396-1653	1
		Neon Pink*	E2396-1654	1
		Neon Yellow*	E2396-1656	1
		Mixed Neon**	E2396-9659	5
	Racks A & B (Complete Workstation) Rack A (without lid) can be inserted into the designated area of Rack B Supplied with lid for Rack B only	Neon Blue*	E2396-9691	1+1
		Neon Green*	E2396-9692	1+1
		Neon Orange*	E2396-9693	1+1
		Neon Pink*	E2396-9694	1+1
		Neon Yellow*	E2396-9696	1+1
		Mixed Neon***	E2396-9669	5+5

* Pack contains one rack with lid. **Pack contains one rack of each colour and five lids. *** Pack contains five Rack A (one of each colour), five Rack B (one of each colour) and five lids for Rack B. Images for illustrative purposes only. Tubes not supplied.

PCR Lock-Rack				
	Description	Colour	Cat. No	Pack Size
	PCR Lock-Rack PP racks which lock together to make larger/smaller racks. Transparent lid fits over two sections locked together. Each section holds 24 x 1.5/2.0 ml tubes, 16 x 0.5 ml tubes and 60 x 0.2 ml tubes. Two sections locked together can hold a 96-Well PCR plate. Dimensions (WxDxH): 133 x 107 x 56 mm (2 sections with lid)	Mixed (contains one section each of blue, green, red, violet, yellow)	E2344-2099	4 Racks, 2 Lids

Images for illustrative purposes only. Tubes not supplied.

Double-Sided PCR Work Rack				
	Description	Colour	Cat. No	Pack Size
	Double-Sided PCR Work Rack Double-sided polypropylene rack with transparent (detachable), hinged lid which has: 168 x 0.2 ml wells AND 12 x 1.5 ml wells on one side, and 40 x 0.5 ml wells AND 12 x 1.5 ml wells on the other Dimensions (WxDxH): 206 x 131 x 54 mm	Mixed (contains one each of blue, green, orange pink, yellow)	E2396-0599	5

Images for illustrative purposes only. Tubes not supplied.