

MiniFuge



Get in a spin!

Small enough to fit on your benchtop for your own personal use. The strip rotor accepts either 16 single or two 8-strip 0.2 ml tubes; the standard rotor accepts six 0.5 ml or 1.5/2.0 ml microcentrifuge tubes. The rotors can be changed easily and quickly.

Features

- Supplied with two, interchangeable rotors for single and strip tubes
- Ideal for quick spin downs and microfiltration
- Starts and stops in seconds



Specifications	
Speed:	6000 rpm
RCF Standard Rotor:	2,000 x g
RCF Strip Rotor:	1,177 x g
Power:	115 V, 60 Hz, 130 mA, or 240 V, 50/60 Hz, 65 mA
Fuse:	0.5 A, 250 V (F)
Dimensions:	150 x 150 x 117 mm

Microcentrifuge		
Cat. No.	Description	Pack Size
N2631-0007	MiniFuge	1

The microcentrifuge is supplied with: One power cord, two rotors ('Standard' for use with 1.5/2.0 ml tubes, and 'Strip' for use with 0.2 ml tubes), six adaptors (to convert the Standard rotor for use with 0.5 ml tubes), and a wrench for changing the rotors.

Images for illustrative purposes only. Tubes not supplied.

*For detailed information please see user manual on our website.